

Folding guide

DryMax Soft a superabsorbent dressing for exuding wounds



In-between toes / fingers




1 Fold 5x5 cm dressing in half.



2 Place the dressing in between toes or fingers.

Heel / elbow



1 Fold 11x20 cm dressing in half widthways.



2 Fold two top corners inwards. Place tape to hold folded corners in place.



3 Open out folded dressing.



4 Place dressing on heel or elbow.

Leg



- 1** Fold 37x56 cm dressing from both sides at the malleolus to cover the foot.



- 2** Continue to fold the dressing around the entire leg.



- 3** Make sure the wound area is covered. Use tape to hold the folded boot shape in place.



- 4** Pass a tube gauze as a sock over the dressing.



- 5** The dressing is in place.

Arm



- 1** Fold 37x56 cm dressing around the arm so that the wound area is covered.



- 2** Thread on a tube gauze over the arm.



- 3** The dressing is in place.

Chest



- 1** Place the 37x56 cm dressing over the chest so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Abdomen



- 1** Place the 37x56 cm dressing over the abdomen so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Back



- 1** Place the 37x56 cm dressing over the back so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.

- 2** Make one cut on each side of a tube gauze to make holes for the arms. Thread on the tube gauze to hold the dressing in place. The dressing is in place.
***Only cut in the gauze tube, not in the dressings.**

Groin, thigh and hip



1

Place the 37x56 cm dressing over the thigh or hip so that the wound area is covered. If necessary, the dressing can be folded and adjusted to the contours of the body.



2

Make one cut of a tube gauze to make hole for the leg. Thread on the tube gauze to hold the dressing in place. The dressing is in place.

***Only cut in the gauze tube, not in the dressings.**

Ordering information

Item	Product size cm	REF	Pcs/box	Absorption capacity*
DryMax Soft	5x5	F60125/10	10	20 ml
DryMax Soft	7.5x7.5	F60011/10	10	60 ml
DryMax Soft	10x10	F60013/10	10	130 ml
DryMax Soft	10x15	F60012/10	10	230 ml
DryMax Soft	12.5x12.5	F60123/10	10	220 ml
DryMax Soft	10x20	F60015/10	10	330 ml
DryMax Soft	15x20	F60014/10	10	400 ml
DryMax Soft	20x20	F60016/10	10	600 ml
DryMax Soft	20x25	F60018/10	10	780 ml
DryMax Soft	20x30	F60017/10	10	960 ml
DryMax Soft	20x40	F60124/10	10	1300 ml
DryMax Soft	37x56	F60019/10	10	3000 ml

WE'RE HERE TO HELP

Email: info@absorbest.com

Phone: +46 (0)494 77 99 50

WWW: absorbest.com



THIS IS DRYMAX SOFT!

Scan the code and find out more

