


DryMax Floor mat XL

The absorbent mat with
anti-slip backing



DryMax Floor mat XL roll



DryMax Floor mat XL mat

DryMax Floor mat XL quickly and efficiently absorbs water-based fluids such as irrigation fluid, blood, urine and other fluids present in examination and operating rooms.

DryMax Floor mat XL is suitable for standing on and can also be placed under a stool.

DryMax Floor mat XL-roll offer flexibility. The XL-roll or any of the mats can be cut into suitable length and form.

Placed on a dry floor the blue anti-slip barrier layer keeps the mat in place and prevents fluid from going through. It prevents slips and falls due to wet floor. Ideal to use under sink area. The mats can, depending on what they have absorbed, be disposed of in the usual manner.

Incorporating DryMax Floor mat XL mats into work routine reduces the time and money spent on cleaning as well as allowing for a faster turnover time.

In the case of heavy liquid flows DryMax Floor mat XL can successfully be used in combination with DryMax Floor mat 2.4 which is a superabsorbent mat, providing exceedingly high absorbency and retention capacity.

Basic 💧 2.0 L/m ²		
REF	Size	Pcs/box
F22201-50	50 cm x 75 cm	55
F22201-100	100 cm x 75 cm	10 x 5 = 50
F22201-160	160 cm x 75 cm	20
F22201-200	200 cm x 75 cm	21
F22201-R	50 m x 100 cm	1
HXL201709	Dispenser for XL roll	1

Medium 💧💧 4.0 L/m ²		
REF	Size	Pcs/box
F22203-100	100 cm x 75 cm	20
F22203-160	160 cm x 75 cm	10
F22203-200	200 cm x 75 cm	11

High 💧💧💧 5.3 L/m ²		
REF	Size	Pcs/box
F22204-100	100 cm x 75 cm	18
F22204-160	160 cm x 75 cm	10
F22204-200	200 cm x 75 cm	9



DISPOSAL
Treat as clinical waste, combustion forms carbon dioxide and water

OPERATION **DryMax**®

Floor mat XL

Safety & Hygiene Control Made Easy

Caution! Place DryMax Floor mat XL on dry floor.



- Examination room
- Operating room
- In sink area/scrub station
- Central Sterile
- Efficient absorption
- Prevents slips and falls
- Faster turnover time
- Reduces time and money spent on cleaning

ABSORBEST

Manufactured by Absorbest AB Sweden
www.absorbest.se | info@absorbest.se