

DryMax Sterile

DryMax Sterile superabsorbent pads lock in large amounts of fluids without dripping. The pads can be used on examination and operating tables to absorb fluids such as saline, urine, amniotic fluids and blood. DryMax Sterile can be easily disposed of and replaced once saturated.


These versatile pads are dual-sided, absorb from both sides and can be used flat, folded or rolled depending on the need of the surgical team.

DryMax Sterile come in two sizes and with a choice of tape along one long side, which enables the pad to be fastened to another surface including the operating drape.

Proper fluid management promotes Health and Safety as it reduces the risk of slips and trips and provides a cleaner and safer working environment.

As DryMax Sterile locks in the contaminated fluids without dripping, it helps reduce the risk of spreading contaminants. In this way, the pad contributes to infection control and patient safety.

Cleaning times and costs are also significantly reduced as fluids are absorbed and retained by the pad.



DISPOSAL
Treat as clinical waste, combustion forms carbon dioxide and water



AREA OF USE

- Arthroscopy
- Knee & shoulder surgery
- Caesarean section
- Amputation
- Abdominal surgery
- Craniotomy
- Burn center

PRODUCT BENEFITS

- Patient & staff safety
- Infection control
- Easy to use
- Faster turnover
- Drier & cleaner environment
- High fluid retention
- No dripping

Item	Product size (cm)	REF	Pcs/box	Packaging (pcs)	Absorption capacity (l)	
					Tap water (l)	NaCl 0.9% (l)
DryMax Sterile Midi	56 x 37	F50161	100	50 x 2 = 100	4 l	1.6
DryMax Sterile Midi Tape	56 x 37	F50163	100	50 x 2 = 100	4 l	1.6
DryMax Maxi	72 x 37	F50241	80	40 x 2 = 80	6 l	2.4
DryMax Maxi Tape	72 x 37	F50243	80	40 x 2 = 80	6 l	2.4

Absorb water-based liquids such as saline, urine, amniotic fluids and blood. In the case of liquids containing bodily fluids and proteins, the absorption capacity will decrease compared to the saline data provided.



DryMax® technology